



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

UREA 6M SOLUTION TESTING METHODS

TABLE OF CONTENTS

1. PURPOSE:.....	3
2. SCOPE:.....	3
3. RESPONSIBILITIES:	3
4. REFERENCES:.....	3
5. EQUIPMENT:.....	3
6. REAGENTS:.....	3
7. ANALYTICAL PROCEDURES:	4

1. PURPOSE:

- 1.1. To provide the Laboratory personnel with a procedure for analyzing Urea 6M In-Process, Finished Goods, and Stability.

2. SCOPE:

- 2.1. Applies to the analysis of Urea 6M In-Process, Finished Goods, and Stability in the Laboratory at all BioSpectra facilities.

3. RESPONSIBILITIES:

- 3.1. The Laboratory Manager, or qualified designee, is responsible for training, maintenance and implementation of this procedure.
- 3.2. The Laboratory Technicians are responsible for compliance with the terms of this procedure. This includes notifying the appropriate personnel if any analyses fail to meet their respective specifications.

4. REFERENCES:

- 4.1. BSI-ATM-0073, Analytical Method of Analysis: Guanidine Thiocyanate, MOPS, and Urea via ICP-MS.
- 4.2. BSI-ATM-0081, Urea 6M Molarity Testing Method via UPLC
- 4.3. BSI-SOP-0098, Balance SOP.
- 4.4. BSI-SOP-0126, Laboratory Notebooks
- 4.5. BSI-SOP-0254, Spectrum Two UATR.
- 4.6. BSI-SOP-0255, XL200 pH/mV/Conductivity Meter SOP.
- 4.7. BSI-SOP-0303, NexION 350X ICP-MS SOP..
- 4.8. BSI-SOP-0345, Endosafe Nexgen-PTS Endotoxin Reader SOP.
- 4.9. BSI-SOP-0350, Anton Paar DMA 35 Portable Density Meter Operation and Calibration SOP.
- 4.10. *ACS, Reagent Chemicals*, current edition.
- 4.11. *Current USP*
- 4.12. *Current EP*

5. EQUIPMENT:

- 5.1. Analytical Balance
- 5.2. Perkin Elmer NexION 350X ICP-MS
- 5.3. Perkin Elmer Spectrum Two UATR
- 5.4. XL200 pH/mV/Conductivity Meter or equivalent
- 5.5. Anton Paar DMA 35 Portable Density Meter
- 5.6. HPLC/UPLC with 195nm Detection Capability

6. REAGENTS:

- 6.1. **LAL Reagent Water:** Purchased Commercially.
- 6.2. **Urea 6M IR Reference Standard:** Utilizing a certified reference standard, prepare by weighing 0.36 grams of reference standard and diluting to 1 mL with purified water (solution may be scaled as needed). Mix to fully dissolve. Prepare at the time of use.

7. ANALYTICAL PROCEDURES:**IN PROCESS TESTING****7.1. MOLARITY** _____ :

- 7.1.1. Refer to the Urea 6M Molarity Testing Method via UPLC, DCN BSI-ATM-0081, for sample preparation, instrumental method parameters and instrument operation for Molarity determination.
- 7.1.2. The weight % will determined via the UPLC analysis, molarity is calculated by the following calculation:

$$7.1.2.1. \text{ Molarity} = \frac{WPs (g)}{100g} \times \frac{Density (g)}{mL} \times \frac{1000mL}{1L} \times \frac{mol}{60.056g}$$

- 7.1.2.1.1. WP_s = Calculated weight percent of the Sample
- 7.1.2.1.2. Density = Measured of the sample
- 7.1.3. **NOTE:** When weighing sample for molarity analysis, perform density analysis at the same time to ensure that temperature is not an impacting factor to the sample being analyzed.

7.2. DENSITY _____ :

- 7.2.1. Analyze the density of the as-is sample utilizing the Anton Paar DMA 35 Portable Density Meter.
- 7.2.2. Perform a water check prior to analysis.
- 7.2.3. Follow the Anton Paar DMA 35 Portable Density Meter Operation and Calibration SOP for instrument operation and calibration, BSI-SOP-0350
- 7.2.4. **NOTE:** When analyzing sample for density, ensure that sample is also being weighed for molarity at the same time to ensure that temperature is not an impacting factor to the sample being analyzed.

FINISHED GOOD TESTING**7.3. APPEARANCE** _____ :

- 7.3.1. Add 50 mL of sample into a Nessler Color Comparison Tube.
- 7.3.2. Add 50 mL of *Purified Water* into a second Nessler Color Comparison Tube.
- 7.3.3. Compare the colors in diffused daylight, viewing vertically against a white background.
- 7.3.4. In order for the sample solution to be colorless, it must have the appearance of *Purified Water*.

7.4. ENDOTOXIN _____ :

- 7.4.1. Pipette 0.140mL of 6M Urea sample into a sterile vessel.
- 7.4.1.1. Note: Total dissolved urea equivalent to 50mg.
- 7.4.1.2. Calculation: 0.140mL*360mg urea /mL = 50mg Urea
- 7.4.2. Dilute sample aliquot to 10mL with LAL reagent water and mix thoroughly.
- 7.4.3. Refer to Endosafe nexgen-PTS Endotoxin Reader SOP, BSI-SOP-0345 for instrument operating instructions.
- 7.4.4. Input a dilution factor of 71.4 into the instrument when prompted.

7.5. IDENTIFICATION (IR) _____ :

- 7.5.1. Prepare a Urea 6M reference standard from a certified reference standard at the time of use.
- 7.5.1.1. Urea 6M = 360mg/mL (360g/L)

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.5.2. Sample will be analyzed as-is
- 7.5.3. Refer to Spectrum Two UATR SOP, BSI-SOP-0254 for instrument operation and standard and sample analysis and comparison.

7.6. **TAMC/TYMC USP <61>** _____ :

- 7.6.1. Aseptically decant and package no less than 50mL of sample into a sterile container and send to MPL laboratories for analysis. The analysis request form should include TAMC/TYMC USP <61> and designate which suitability method to be used on the Analysis Request Form (ARF). A new suitability should be requested if a reference suitability is unavailable.

7.7. **MOLARITY** _____ :

- 7.7.1. Refer to the Urea 6M Molarity Testing Method via UPLC, DCN BSI-ATM-0081, for sample preparation, instrumental method parameters and instrument operation for Molarity determination.
- 7.7.2. The weight % will determined via the UPLC analysis, molarity is calculated by the following calculation:

$$7.7.2.1. \text{ Molarity} = \frac{WPs (g)}{100g} \times \frac{\text{Density} (g)}{mL} \times \frac{1000mL}{1L} \times \frac{mol}{60.056g}$$

7.7.2.1.1. WP_s = Calculated weight percent of the Sample

7.7.2.1.2. Density = Measured of the sample

- 7.7.3. **NOTE:** When weighing sample for molarity analysis, perform density analysis at the same time to ensure that temperature is not an impacting factor to the sample being analyzed.

7.8. **DENSITY** _____ :

- 7.8.1. Analyze the density of the as-is sample utilizing the Anton Paar DMA 35 Portable Density Meter.
- 7.8.2. Perform a water check prior to analysis.
- 7.8.3. Follow the Anton Paar DMA 35 Portable Density Meter Operation and Calibration SOP for instrument operation and calibration, BSI-SOP-0350.
- 7.8.4. **NOTE:** When analyzing sample for density, ensure that sample is also being weighed for molarity at the same time to ensure that temperature is not an impacting factor to the sample being analyzed.

7.9. **pH @ 25°C** _____ :

- 7.9.1. Calibrate the XL200 pH/mV/Conductivity meter with pH 4, 7, and 10 buffer standards.
- 7.9.2. Measure the neat sample and record the results at $25 \pm 2^\circ\text{C}$.
- 7.9.3. Refer to XL200 pH/mV/Conductivity SOP for instrument operation and sample analysis, BSI-SOP-0255.

7.10. **TRACE ELEMENTS** _____ :

- 7.10.1. Refer to Analytical Method of Analysis: Guanidine Thiocyanate, Mops, and Urea via ICP-MS (BSI-ATM-0073) for sample preparation and analysis.