

2/16/2026

Version 1.2

55 GAL/200L BLACK POLY DRUM SPECIFICATIONS



55-Gallon Black Poly Drum	
Test	Specifications
Appearance and Color	Black, closed head drum with white click lock caps
Outside Height (Nominal)	34.8"
Outside Diameter, Top (Nominal)	23.3"
Weight (Nominal)	9.52 kg. (including caps)
Capacity	55.0 U.S. Gallons (USG)
UN Rating	1H1/Y1.9/150/YR

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.