




**Product Data**  
 Life Latch® Drum  
 12.0 Gallon - Tamper-proof

**Cover Specification**

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Weight	1.85 lbs. / 840g
Height	2 1/2"
Min Wall Thickness	0.100"
Bottom Outside Diameter	16 3/8"
Gasket	Neoprene / Available w/o Gasket
Fitments	Tamper-proof Zip Tie Window <i>Please inquire about custom fittings</i>
Standard Color	White <i>Please inquire about custom colors</i>
Shipping Quantity	180 per pallet



**Pail Specification**

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Handles	NA
Weight	5.91 lbs. / 2680 g
Exterior Height	21 1/4" / With Cover 21.35"
Min. Wall Thickness	0.14"
Inside Diameter - Top	14 1/8"
Outside Diameter - Bottom	13"
Standard Color	White <i>Please inquire about custom colors</i>
Shipping Quantity	54 Units per Pallet 2,160 Units per Truckload
Print Area	12" x 41 1/4"
Top Load Rating	450 lbs.

