

2/16/2026

Version 1.2

4-UNIT DESICCANT SPECIFICATIONS



4-Unit Tyvek Desiccant	
Test	Specifications
Desiccant Name	4-Unit Tyvek Desiccant
Desiccant Type	Desi Pak®
Appearance and Color	Sealed White Tyvek® package with Black or Blue print
Stamp Markings	Specification MIL-D-3464 Type I & II
Seal Integrity	Physically inspect the Seal and ensure it is fully intact with the seal firmly applied

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.