

2/16/2026

Version 1.2

## 1KG BOTTLE SPECIFICATIONS



1kg HDPE Lab Screw-Top Bottle	
Test	Specifications
Appearance and Color	White plastic bottle with wide mouth and recessed label panel
Lid	Blue, Tamper-Evident Screw Cap
Capacity	50 oz / 1500 mL
Outside Height	8.25 inches $\pm$ 0.5 inches including lid
Outside Width	4.25 inches $\pm$ 0.5 inches
UN Rating	3H1/X4/S

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director  
Business Integration & Retention  
BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*