

500G BOTTLE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0033



500g HDPE Lab Screw-Top Bottle	
Test	Specifications
Appearance and Color	White plastic bottle with wide mouth and recessed label panel
Lid	Blue, Tamper-Evident Screw Cap
Capacity	33 oz / 1000 mL
Outside Height	8.75 inches \pm 0.5 inches including lid
Outside Width	3.25 inches \pm 0.5 inches
UN Rating	3H1/X2.5/S

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.