

# BIOBUFFER SOLUTIONS

100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

December 16<sup>th</sup>, 2024

Revision 1

## RESIDUAL SOLVENTS STATEMENT

### L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Qualified Grade is manufactured using the solvent Ethyl Acetate, and is stated to comply with the ICH Q3C (R8) Guideline for Residual Solvents required concentration limit of 5000ppm for a Class 3 residual solvent. BioSpectra has evaluated the material supply through the Supplier Qualification Program.

Current Product Number
LARH-6201

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)



**Cassie Baun**

Senior Compliance Specialist