

2/16/2026

Version 1.4

4x1 GALLON (3.78L) BOTTLE SPECIFICATIONS



| 4x1 Gallon Industrial Shippers | |
|--------------------------------|--|
| Test | Specifications |
| Appearance and Color | 4x1 Gallon Non-Fluorinated HDPE Natural Round Bottle (Includes carton, corrugated pads, poly bags, locking ties and divider) |
| Lid | Polypropylene Cap with Polyethylene Foam Liner |
| Capacity | 4 Bottles x 1 Gallon Each |
| Outside Height | 11.64 inches \pm 0.5 inches |
| Outside Width | 6.10 inches \pm 0.5 inches |
| UN Rating | 4G/Y26/S |

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.