


UN Certification: 1H2/Y30/S & 1H2/X11.5/S


Cover Specification | Part No. 600137

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding
Weight	1.03 lbs / 470g
Height	2.50"
Min Wall Thickness	.090"
Bottom OD	12.69"
Gasket	Neoprene / Available w/o Gasket
Standard Color	White, Red 186, Blue 300 <i>Please inquire about custom colors</i>
Case Quantity	20 per case / 240 per pallet


Pail Specification | Part No. 500071

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding / Plastic Handle
Weight	2.36 lbs / 1071 g
Min. Wall Thickness	.090"
Exterior Height	15.25" With Cover 16.04"
Top Inside Diameter	11.25"
Bottom Outside Diameter	10.27"
Shipping Quantity	4 Pails and Covers Per Case 120 Units per Pallet 4,320 Units per Truckload
Standard Color	White <i>Please inquire about custom colors</i>
Print Area	9.5" x 33"
Heat Transfer Area	8.5" x 31"
Top Load Rating	425 lbs


 RECEIVED
 APR 28 2020
