

10 Liter opti rib® can HDPE, 450 g, DIN 50, UN

Art.-Nr.	722 006 30	722 006 12	722 006 32
Ident-Nr.:	2000095642	2000095902	2000095811
Colour:	natural	white	blue

Description

with blow-moulded Handle
 non-residue emptying
 label field on both sides
 stackable
 with visibility strip (optional)

Dimension:

Length 230 mm ±3mm
 Width 190 mm ±3mm
 Height 307 mm ±4mm
 Weight 450 g + 4 g
 Overflow capacity 10,8 Liter
 Neck Opening Ø 40,5 mm -0,5mm

Label field:

Height 168,0 mm
 Width 160,0 mm

Imprints:	Front	Back	Shoulder	Bottom
	Scale 1 - 9 l	Braille imprint	UN certification	Kautex
	Braille imprint		PE-Recycling symbol	10 Liter
			Production date	Mould No.

Packing units	10 pieces / package	2 package / layer	4 layers / pallet
Amount on pallet	80 pieces	180 pieces	Euro pallet (1200 x 800 mm)
Height of pallet	1336 mm	2796 mm	

Approval: 3H1/X1.9/250/ .. /D/ BAM 14715-KWA

This can is certified for the filling of dangerous materials according to the European Dangerous Goods Regulations. Certification applies to the following substances:

Model substance	Density
Water	1,9 g /cm ³
Wetting Agents	1,2 g /cm ³
Acetic Acid	1,4 g / cm ³
n-Butyl Acetate	1,2 g / cm ³
White spirit	1,2 g / cm ³
Nitric Acid 55 %	1,2 g /cm ³

The UN-certification is only valid with approved closures. (see details in the certificate of approval)

approved closures	drawing no.	raw material	supplier
SKGB 45261 R04	01-122/3k	L	GFV Verschlusstechnik GmbH & Co. KG
SKGS 45271 R05	93-303/1k	L	GFV Verschlusstechnik GmbH & Co. KG
SKUB 45261 R04	04-117/1k	L	GFV Verschlusstechnik GmbH & Co. KG
SKUB 45261 R04 EC	05-107/1k	L	GFV Verschlusstechnik GmbH & Co. KG
SK 45/26 MAB MDR	ENG-00-000174 g	I	BERICAP GmbH & Co. KG
SK 45/26 MAB MDR MPV D15	ENG-00-016609 d	I	BERICAP GmbH & Co. KG
SK 45/26 MAB MDR GE 4	ENG-00-004665 e	I	BERICAP GmbH & Co. KG
SK 45/26 MAB MDR MPV D17	ENG-00-016585 g	I	BERICAP GmbH & Co. KG
SK 55/27 MAB MDR (DIN)	ENG-00-008832 g	L	
SK 55/27 MAB MDR MPV D15 (DIN)	ENG-00-014431 d	L	

Product Information



820	P050-010 Zust. 1	L	
820-VE	P050-011 Zust. 1	L	
SV 51 OV-A PE	SV51 OV-A PE	L	
SV 51 OV-A VE PE	SV51 OV-A VE PE	L	
SK 51/16 OV TC flach	V 5134	O	

No responsibility is taken for the correctness of this product information. This document is subject to change without prior written notice.