

2/16/2026

Version 1.2

100G BOTTLE SPECIFICATIONS



100g HDPE Lab Screw-Top Bottle	
Test	Specifications
Appearance and Color	White plastic bottle with wide mouth and recessed label panel
Lid	Blue, Tamper-Evident Screw Cap
Capacity	8 oz / 250 mL
Outside Height	5.25 inches \pm 0.5 inches including lid
Outside Width	2.25 inches \pm 0.5 inches
UN Rating	3H1/X0.8/S

For further information, please contact marketing@biospectra.us.

Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.