

## 50KG POLY DRUM SPECIFICATIONS



Blue 30 Gallon / 50kg Poly Drum		
Item	Test	Specifications
Drum	Appearance and Color	Light Blue Plastic Drum
	Nominal Capacity	30.0 Gallons 113.6 Liters Drum Stamp: 114 Liters
	Height (Not including lid)	Body Height: 29.8"
	Height (Including lid)	Overall Height: 30.1"
	OD (Top)	Nominal OD: 21.5"
	OD (Bottom)	Bottom OD: 17.2"
	ID (Top)	Opening ID: 19.2"
	Hand Grip	Two (2) shallow hand grips located on base of drum
	UN Markings on Certificate	1H2/Y180/S/YR at minimum
Lid	Gasket	Crème Colored Gasket
	Gasket Type	612
	Gasket Height	0.345-0.385"
	Rim Diameter	20.77" ± 0.03"
	Skirt Diameter	Skirt Diameter: 19.22" ± 0.03"
	Appearance and Color	Light Blue with Contours and Recycling Symbol 2 HDPE
Lock Band	Description	20" Open Head Loc Band – Natural – 69.5" Dimension – UV – No Ribs Open head handle – Natural – UV – Standard
	Color	Natural

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director of Commercial Operations**  
**BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*