

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	8-Aug-2022	8-Aug-2025	: Date of Next Review
Prepared By:	Amy Hosein	BSI-COA-0001 v.6.0	: Supersedes
QA/QC Approval:	Brett Brink	Carissa McCollian	: Management Approval
Reason for Revision:	See Revision History in MasterControl.		

CERTIFICATE OF ANALYSIS

TRIS HCl

BIO EXCIPIENT GRADE / NEW CODE THCL-3220-10

(HISTORICAL CODE TH3220-K010)

LOT: THCL-0123-00264

 $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3 \cdot \text{HCl}$ * F.W. 157.60 g/mol. * CAS# 1185-53-1

Manufacturing Date: 08/30/23 Expiration Date: 08/31/26

Manufacturing Site: 1474 Rockdale Lane, Stroudsburg, PA 18360

Packaging Date: 10/09/23 Packaging Site: 100 Majestic Way, Bangor PA, 18013

ANALYSIS	SPECIFICATION	TEST RESULT	
Absorbance (1M)	280 nm	0.06 a.u. max.	0.01 a.u.
Appearance and Color	White / Crystals		Passes Test
Assay	99.5% min.		99.9%
Enzymes	DNase	None Detected	None Detected
	RNase	None Detected	None Detected
	Protease	None Detected	None Detected
Heavy Metals	2 ppm max.		< 2 ppm
Identification (IR)	Passes Test		Passes Test
Insoluble Matter	0.001% max.		0.001%
Karl Fischer	0.5% max.		0.4%
Melting Range	150 - 153 °C		150 - 151 °C
pH (0.5M)	4.0 - 5.0		4.3 @ 23.2 °C
pK _a	8.0 - 8.4		8.2
Residue on Ignition	0.1% max.		< 0.1%
Solubility 35%	Passes Test		Passes Test
Trace Metals	Arsenic (As)	1 ppm max.	< 0.45 ppm
	Calcium (Ca)	1 ppm max.	< 0.60 ppm
	Copper (Cu)	1 ppm max.	< 0.15 ppm
	Iron (Fe)	1 ppm max.	< 0.30 ppm
	Lead (Pb)	1 ppm max.	< 0.30 ppm
	Magnesium (Mg)	1 ppm max.	< 0.60 ppm

COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: BSI-ATM-0002

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

RESIDUAL SOLVENTS: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by: M. Shapiro Date: 10/17/23 Job Title: QA Mater. Disp. Tech. III

Reviewed by: Carrie Albert Date: 10/17/23 Job Title: Assoc. Director of Quality