

TRIS

MEETS USP, EP, JPC, ChP

CAS #: 77-86-1

Formula: $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3$

F.W.: 121.14 g/mol

TRIS-6202

BIO QUALIFIED GRADE

ANALYSIS	SPECIFICATIONS
Appearance of Solution (EP)	Passes Test
Assay, Dried Basis (USP/EP/ChP/JPC)	99.0 – 100.5%
Chloride (ChP/EP)	< = 100 ppm
Clarity and Color of Solution (ChP/JPC)	Passes Test
Heavy Metals (Pb) (ChP/JPC)	< = 8 ppm
Identification, IR (USP-A/ChP-3/EP-C)	Conforms to Reference Standard
Identification A (EP), Alkalinity (ChP), pH (5%) (USP/EP)	10.0 – 11.5
Identification A (JPC)	Passes Test
Identification B (USP/ChP-1)	Passes Test
Identification B (EP), Melting Range (USP/ChP/JPC)	168 – 172 °C
Identification B (JPC)	Passes Test
Identification C (USP)	Passes Test
Identification D (EP/ChP-2)	Conforms to Reference Standard
Loss on Drying (USP/ChP/EP/JPC)	< = 0.5%
pH (1 in 100) (JPC)	10.3 – 10.7
Related Substances (ChP/EP)	< = 1.0%
Residue on Ignition/Sulfated Ash (USP/ChP/EP/JPC)	< = 0.1%
Trace Metals	Arsenic (As) (JPC) <= 1.6 ppm
	Iron (Fe) (ChP/EP) < = 10 ppm
	Nickel (Ni) (ChP) < = 15 ppm

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

General Product Overview

Tromethamine (Tris) is commonly used as a buffer and excipient in downstream biological applications during the manufacture of drug products. Tris is considered a Good's buffer because it has low UV absorptivity, minimal reactivity, stable pH and is soluble in water.

Industry Application

Suitable for use in biological and biotech chemical process applications from R&D through scale cGMP production.

Key Product Features

- Appears as a white crystalline powder
- Manufactured in accordance with the ISO Certified Quality Management System
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: Trometamol; THAM; Tris (Hydroxyethyl) Aminomethane; Trisaminol; Trismethylolaminomethane; 2-Amino-2(Hydroxymethyl) propane-1,2-diol; Tromethamine

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Please inquire for information regarding shelf life.

Package Sizes

1kg, 5kg, 10kg, 25kg, 50kg

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.