

12/18/2025

Version 1.3

6X1L BOTTLE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0030



6x1 Liter Non-Fluorinated HDPE Bottle Combo Pack	
Test	Specifications
Appearance and Color	Natural HDPE Bottles in UN Kraft Carton – Non-Fluorinated
Capacity	6 Bottles x 1 Liter Each
Outside Height	10.455 inches – 10.615 inches
Air Insert (Bottom of bottle)	Arrow pointed at A
UN Rating	4G/Y13.9/S / 4G/Y14.1/S

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.