



100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

## HDPE CARBOY COMPONENT SPECIFICATIONS

2.6 Gallon HDPE Carboy	
Test	Specifications
Capacity	2.6 Gallons
Material	Natural HDPE Jerry Can
Appearance and Color	White, rectangular with built-in handle
Outside Length (Internal Tolerance)	9.1 inches $\pm$ 1 inch
Outside Width (Internal Tolerance)	7.5 inches $\pm$ 1 inch
Outside Height (Internal Tolerance)	12.1 inches $\pm$ 1 inch
Closure	Continuous Thread Round Screw Top

Supplier Part #: 1248B07-B, or 1248B07

Manufacturer Part #: 2000095641 or 2000095642 (Identical Item Per Supplier)

Supplied by: Berlin Packaging, LLC

Manufactured by: Kautex

2.6 Gallon HDPE Carboy Cap	
Test	Specifications
Material	50mm Black HDPE Cap
Appearance and Color	Black Round Tamper Evident Screw Cap

Supplier Part #: 235649, or 1248C03

Manufacturer Part #: 2000095946

Supplied by: Berlin Packaging, LLC

Manufactured by: Kautex

<b>5 Gallon HDPE Carboy</b>	
<b>Test</b>	<b>Specifications</b>
Capacity	5 Gallons
Material	Natural HDPE
Appearance and Color	White, rectangular with built-in handle, vent lid and spigot adapter
Outside Length (Internal Tolerance)	10.00 inches $\pm$ 0.5 inch
Outside Width (Internal Tolerance)	9.75 inches $\pm$ 0.5 inch
Outside Height (Internal Tolerance)	15 inches $\pm$ 0.5 inch including lid
Mouth Opening (Internal Tolerance)	2 and 3/8 inches $\pm$ 0.5 inches

Part #: H002/H003

Supplied by: Container &amp; Packaging Supply

Manufactured by: Zacros America Inc. Hedwin Div./ Streamline Plastics

<b>20 Liter (5.3 Gallon) Nalgene HDPE Carboy</b>	
<b>Test</b>	<b>Specifications</b>
Capacity	20 Liters (5.3 Gallons)
Material	HDPE Carboy, 100/415 PP Cap
Appearance and Color	White, rectangular with metal handle and Round White Screw Cap
Outside Length (Internal Tolerance)	12.5 inches $\pm$ 0.5 inch
Outside Width (Internal Tolerance)	8.75 inches $\pm$ 0.5 inch
Outside Height (Internal Tolerance)	15.50 inches $\pm$ 0.5 inch
Closure	Continuous Thread Round Screw Top

Catalog #: 02-961-10B

Manufacturer Part #: 2211-0050

Supplied by: Fisher Scientific

Manufactured by: Thermo Scientific