

2/16/2026

Version 1.2

TYVEK LINER SPECIFICATIONS



Tyvek Liner Component Specifications	
Test	Specifications
Appearance and Color	White with no punctures
Width	36.0" \pm 0.5"
Length (Internal Tolerance)	60.0" \pm 1.0"

For further information, please contact marketing@biospectra.us.

Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.