

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	01-Mar-2021	01-Mar-2024	: Date of Next Review
Prepared By:	Jared L Lobb	16-001173 v.4.1	: Supersedes
QA/QC Approval:	Jaron Hughes	Wendy Santay	: Management Approval
Reason for Revision:	See Revision History in ensur.		

CERTIFICATE OF ANALYSIS

MES MONOHYDRATE

BIO EXCIPIENT GRADE / NEW CODE MESM-3220-93

(HISTORICAL CODE ME3220-G500)

LOT: MESM-0122-00077

$C_6H_{13}NO_4S \cdot H_2O$ Δ F.W. 213.3 g/mol. Δ CAS# 145224-94-8

Manufacturing Date: 2/24/22 Retest Date: 2/29/24

Manufacturing Site: 100 Majestic Way, Bangor PA, 18013

Packaging Date: 6/7/22 Packaging Site: 100 Majestic Way, Bangor PA, 18013


ANALYSIS	SPECIFICATION	TEST RESULT
Absorbance (1M)	260 nm 280 nm	0.1000 a.u. max. 0.1000 a.u. max.
Appearance and Color	White / Crystals	0.0026 a.u. 0.0016 a.u.
Assay	99.0% min.	White / Crystals
Chloride	0.005% max.	100.0%
Color (1M, Alkaline)	Colorless	< 0.005%
Enzymes	DNase	Colorless
	RNase	None Detected
	Protease	None Detected
Heavy Metals (as Pb)	2 ppm max.	None Detected
Identification (IR)	Passes Test	None Detected
Loss on Drying @ 130°C	7 – 10%	None Detected
pH (5% Soln.)	3.1 – 3.5	9%
pH (0.5M)	2.5 – 4.5	3.5
pK _a	5.9 – 6.3	3.2
Residue on Ignition	0.05% max.	6.1
Solubility (5%)	Passes Test	< 0.01%
Sulfate	0.005% max.	Passes Test
Trace Elements	Arsenic (As)	< 0.005%
	Copper (Cu)	2 ppm max.
	Iron (Fe)	2 ppm max.
	Lead (Pb)	2 ppm max.


COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: 16-001016

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

OVI STATEMENT: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by:  Date: 6/8/22 Job Title: QA specialist

Reviewed by:  Date: 6/8/22 Job Title: QA Manager