

2/16/2026

Version 1.2

STANDARD HEAT TREATED PALLET SPECIFICATIONS



48" x 42" Standard Wooden Pallets, Heat Treated	
Test	Specifications
Appearance and Color	Light Brown
Planks	Brown wood planks, 3 stringers
Length (Nominal)	48.0"
Width (Nominal)	42.0"
Height (Nominal)	4.75"
Stamp	Heat Treated

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.