

GUANIDINE THIOCYANATE

CAS #: 593-84-0

 Formula: C₂H₆N₄S

F.W.: 118.16 g/mol

GTHI-5221
BIO ULTRA GRADE

ANALYSIS		SPECIFICATIONS
Appearance and Color		White Crystals
Assay, Dried Basis		> = 99.0%
Enzymes	RNase DNase	None Detected None Detected
Identification, IR		Conforms to Reference Standard
Solubility (100mg/mL)		Clear, Colorless
Trace Metals	Barium (Ba)	< = 5 ppm
	Iron (Fe)	< = 5 ppm
	Lead (Pb)	< = 5 ppm
	Potassium (K)	< = 50 ppm
	Sodium (Na)	< = 5000 ppm

Key Product Features

- Appears as white crystalline powder
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: Guanidine Isothiocyanate

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Please inquire for information regarding shelf life.

Package Sizes

1kg, 5kg, 10kg, 25kg, 50kg

Industry Application

Suitable for use in biological and biotech chemical process applications from R&D through scale production.

General Product Overview

Guanidine Thiocyanate is a chaotropic agent. As one of the strongest protein denaturants, it inactivates nucleases more efficiently than Guanidine Hydrochloride. Guanidine Thiocyanate is preferred for the purification of RNA because it dissociates the RNA into its nucleic acids and protein forms. Higher molarity solutions of Guanidine Thiocyanate irreversibly inactivate RNase. Guanidine Thiocyanate causes proteins to dissolve readily, thus disintegrating cellular structures.

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.