

# BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

Effective Date:	8-Aug-2022	8-Aug-2025	: Date of Next Review
Prepared By:	Army Hosein	BSI-COA-0001 v.6.0	: Supersedes
QA/QC Approval:	Brett Brink	Carissa McCollan	: Management Approval
Reason for Revision:	See Revision History in MasterControl.		

## CERTIFICATE OF ANALYSIS

### TRIS HCl

### BIO EXCIPIENT GRADE / NEW CODE THCL-3220-10

(HISTORICAL CODE TH3220-K010)

LOT: THCL-0123-00088

$\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3 \cdot \text{HCl}$  \* F.W. 157.60 g/mol. \* CAS# 1185-53-1

Manufacturing Date: 2/13/23      Expiration Date: 2/28/26

Manufacturing Site: 100 Majestic Way, Bangor PA, 18013

Packaging Date: 2/28/23 Packaging Site: 100 Majestic Way, Bangor PA, 18013

ANALYSIS	SPECIFICATION	TEST RESULT
Absorbance (1M)	280 nm 0.06 a.u. max.	0.01 a.u.
Appearance and Color	White / Crystals	Passes Test
Assay	99.5% min.	99.8%
Enzymes	DNase	None Detected
	RNase	None Detected
	Protease	None Detected
Heavy Metals	2 ppm max.	<2 ppm
Identification (IR)	Passes Test	Passes Test
Insoluble Matter	0.001% max.	<0.001%
Karl Fischer	0.5% max.	0.3%
Melting Range	150 – 153 °C	151-152°C
pH (0.5M)	4.0 – 5.0	4.3 @ 23.0°C
pK <sub>a</sub>	8.0 – 8.4	8.2
Residue on Ignition	0.1% max.	<0.1%
Solubility 35%	Passes Test	Passes Test
Trace Metals	Arsenic (As)	1 ppm max. <0.45 ppm
	Calcium (Ca)	1 ppm max. <0.60 ppm
	Copper (Cu)	1 ppm max. <0.15 ppm
	Iron (Fe)	1 ppm max. <0.30 ppm
	Lead (Pb)	1 ppm max. <0.30 ppm
	Magnesium (Mg)	1 ppm max. <0.60 ppm

COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: BSI-ATM-0002

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

RESIDUAL SOLVENTS: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by: *Jan Singh* Date: 3/13/23 Job Title: QA Mater. Disp. Supervisor  
Assoc. Director

Reviewed by: *Camille* Date: 3/13/23 Job Title: of Quality