

Product size	Outer packaging volume	Packaging Weight	Outer packaging description	Inner packaging size	Inner packaging description	Inner packaging seal	Outer Packaging seal
10 - 100 GRAMS	9 x 7 x 4" IN	0.5 kg	CORRUGATED BOX	5 x 8" IN	4 MIL MINIGRIP RECLOSABLE BAGS	8" CABLE TIE	N/A
500 GRAMS	9 x 7 x 4" IN	0.5 kg	CORRUGATED BOX	10 x 12" IN	4 MIL BAGS ON A ROLL	8" CABLE TIE	N/A
1 - 2 KILOGRAMS	10 x 7 x 7" IN	1.00 kg	CORRUGATED BOX	10 x 12" IN	4 MIL BAGS ON A ROLL	8" CABLE TIE	N/A
3 - 20 KILOGRAMS	15 GALLON	2.5 kg	FIBER DRUM	(2 BAGS) 24 x 36" IN	3 MIL BOTTOM SEAL BAG	8" CABLE TIE	N/A
21 - 70 KILOGRAMS	55 GALLON	7.5 kg	FIBER DRUM	(2 BAGS) 38.5 x 72" IN	3 MIL FULL GAUGE BOTTOM SEAL BAG	14" CABLE TIE	14" CABLE TIE
100 MLS	9 x 7 x 4" IN		CORRUGATED BOX	100 MLS	GLASS VIAL	CARDBOARD DIVIDER	N/A
500 MLS	9 x 7 x 4" IN		CORRUGATED BOX	1 LITRE	PLASTIC BOTTLE	N/A	N/A
1 LITRE	9 x 7 x 4" IN		CORRUGATED BOX	1 LITRE	PLASTIC BOTTLE	N/A	N/A
10 LITRES	12 x 12 x 18" IN		MULTI-DEPTH CORRUGATED BOX	10 LITRE	BLUE/YELLOW CUBATAINER	N/A	N/A
20 LITRES	12 x 12 x 18" IN		MULTI-DEPTH CORRUGATED BOX	20 LITRE	BLUE/YELLOW CUBATAINER	N/A	N/A
200 LITRES	200 L	10 kg	BLUE PLASTIC DRUM	N/A	N/A	N/A	N/A
1000 LITRES	275 GALLON	60 kg	TOTE (ECOBULK #4038263)	N/A	N/A	48x40x48 THERMAL PALLET COVER (IN WINTER)	N/A
1000 LITRES	275 GALLON	60 kg	TOTE (RECOBULK #1332007)	N/A	N/A	N/A	N/A

Shipping Weights and Dimensions: (Only 1 type of Skid)

Wooden Skid Weight: 19 kgs

Wooden Skid Weight: 24 lb.

Dimensions: L: 48' x W: 44" x H: 5 "

Dimensions:

SMALL DRUM:

Skid Height + Drum Height = 22 inches

Skid Width = 44 inches.

Skid Length = 48 inches

LARGE DRUM:

Skid Height + Drum Height = 39.8 inches

Skid Width = 44 inches.

Skid Length = 48 inches

CUBATAINER:

Skid Height + Cubtainer = 20.375 inches

Skid Width = 44 inches.

Skid Length = 48 inches