

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

TOTE COMPONENT SPECIFICATIONS

330-Gallon Tote	
Test	Specifications
Outside Width (Internal Tolerance)	40" ± 1.0"
Outside Length (Internal Tolerance)	48" ± 1.0"
Outside Height (Internal Tolerance)	53" ± 1.0"
Fill Opening Screw Cap (Internal Tolerance)	6" ± 0.5"
Discharge Opening (Internal Tolerance)	2" ± 0.5"
Capacity	330 gallons
Pallet Type	Steel Frame Pallet

Part #: 1208233

Supplied by: Schuetz

Manufactured by: Schuetz

275-Gallon Tote	
Test	Specifications
Outside Width (Internal Tolerance)	39" ± 1.0"
Outside Length (Internal Tolerance)	47" ± 1.0"
Outside Height (Internal Tolerance)	46" ± 1.0"
Fill Opening Screw Cap (Internal Tolerance)	6" ± 0.5"
Discharge Opening (Internal Tolerance)	2" ± 0.5"
Capacity	275 gallons
Pallet Type	Steel Frame Pallet

Part #: 939544

Supplied by: Schuetz

Manufactured by: Schuetz

275-Gallon Tote	
Test	Specifications
Outside Width (Internal Tolerance)	39" ± 1.0"
Outside Length (Internal Tolerance)	47" ± 1.0"
Outside Height (Internal Tolerance)	46" ± 1.0"
Fill Opening Screw Cap (Internal Tolerance)	6" ± 0.5"
Discharge Opening (Internal Tolerance)	2" ± 0.5"
Capacity	275 gallons
Pallet Type	Plastic Skid Pallet

Part #: 4020528
 Supplied by: Schuetz
 Manufactured by: Schuetz