

BIOBUFFER SOLUTIONS

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

December 16th, 2024

Revision 1

SUPPLY CHAIN STATEMENT

L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Qualified Grade is manufactured and packaged by BioSpectra's approved supplier in accordance with the approved supplier's ISO 9001:2015 certified management system. Additional packaging of L-Arginine Hydrochloride, Bio Qualified Grade is performed at BioSpectra's Bangor, PA Facility: 100 Majestic Way, Bangor, PA 18013. The storage and shipping of L-Arginine Hydrochloride, Bio Qualified Grade occurs in BioSpectra's Supply Chain Center: 51 North 3rd Street, Stroudsburg, PA 18360.

Current Product Number
LARH-6201

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist