

## QUALITY CONTROL TESTING AND VALIDATION STATEMENT

### L-Histidine Monochloride Monohydrate

L-Histidine Monochloride Monohydrate, Bio Pharma Grade for BioBuffer Solutions Product line is analyzed by BioSpectra's Quality Control laboratory. Analytical methods utilized for the analysis of L-Histidine Monochloride Monohydrate, Bio Pharma Grade for BioBuffer Solutions Product line are validated as applicable in accordance with BioSpectra's Analytical Methods Validation Master Plan. The L-Histidine Monochloride Monohydrate testing methods listed in the table below are validated in accordance with USP General Chapter requirements.

Current Product Number
LHMM-4250

Analysis Name
Assay (different calculations to report dried substance and anhydrous basis)
Endotoxin
Trace Metals to report: Heavy Metals and Iron
Water

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)



**Cassie Baun**

Senior Compliance Specialist