

2/16/2026

Version 1.2

24" x 36" POLY LINER SPECIFICATIONS



24" x 36" Poly Liner	
Test	Specifications
Appearance and Color	Clear with no punctures
Dimensions (Nominal)	24" x 36"
Thickness	0.0027" minimum
Materials of Construction (Film)	Hexene LLDPE Nonslip
FDA Statement	21 CFR 177.1520
UATR	Passes Test; when compared to a previously approved poly liner for this part number

For further information, please contact marketing@biospectra.us.

Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.