

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
BTRI-3250-01	Bis Tris, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
BTRI-3250-05	Bis Tris, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-10	Bis Tris, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-25	Bis Tris, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-3250-92	Bis Tris, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
BTRI-3250-93	Bis Tris, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
BTRI-4250-01	Bis-Tris Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
BTRI-4250-05	Bis-Tris Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-10	Bis-Tris Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-25	Bis-Tris Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-4250-50	Bis-Tris Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
BTRI-4250-92	Bis-Tris Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
BTRI-4250-93	Bis-Tris Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
CSMH-3250-01	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		1.25 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-05	Cysteamine HCl (2-MEA), 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-10	Cysteamine HCl (2-MEA), 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-25	Cysteamine HCl (2-MEA), 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
CSMH-3250-90	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		3.50 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-91	Cysteamine HCl (2-MEA), 50 G	In development	50.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-92	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-94	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3350-10	Cysteamine HCl (2-MEA), 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3350-25	Cysteamine HCl (2-MEA), 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
D010-3201-01	Dextran Powder MW 10,000, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
D010-3201-30	Dextran Powder MW 10,000, 30 KG	Active	30.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXSE-4250-01	Dextran Sulfate, 8000, GMP, 1kg	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-05	Dextran Sulfate, 8000, GMP, 5kg	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-10	Dextran Sulfate, 8000, GMP, 10kg	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-15	Dextran Sulfate, 8000, GMP, 15kg	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-25	Dextran Sulfate, 8000, GMP, 25kg	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-92	Dextran Sulfate, 8000, GMP, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-93	Dextran Sulfate, 8000, GMP, 20g	Active	20.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-94	Dextran Sulfate, 8000, GMP, 500g	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-01	Dextran Sulfate Sodium Salt MW 8000 Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-05	Dextran Sulfate Sodium Salt MW 8000 Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-10	Dextran Sulfate Sodium Salt MW 8000 Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-25	Dextran Sulfate Sodium Salt MW 8000 Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4252-50	Dextran Sulfate Sodium Salt MW 8000 Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXSE-4252-92	Dextran Sulfate Sodium Salt MW 8000 Biotech, 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-93	Dextran Sulfate Sodium Salt MW 8000 Biotech, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-5201-92	Dextran Sulfate, 8000, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
GALP-3224-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3224-25	D-Galactose, Plant Derived, 20kg	In development	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3250-05	D-Galactose, Plant Derived, 5kg	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GALP-3250-10	D-Galactose, Plant Derived, 10kg	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-25	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3250-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3251-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3251-10	D-Galactose, Plant Derived, 10kg	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3251-21	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3251-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3251-50	D-Galactose, Plant Derived, 50kg	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3252-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3252-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3350-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3350-25	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3350-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3351-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-50	D-Galactose, Plant Derived, 50kg	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-05	Galactose (D) Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-06	Galactose (D) Biotech, 1 KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-4250-10	Galactose (D) Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-21	Galactose (D) Biotech, 20 KG	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-4250-27	Galactose (D) Biotech, 25 KG	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-50	Galactose (D) Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-92	Galactose (D) Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GHCL-3101-02	Guanidine HCl 6M Solution, 1L	Active	1.00	L	1.00	L	Natural HDPE bottle			
GHCL-3101-72	Guanidine HCl 6M Solution, 200 L	Active	200.00	L	55.00	GAL	Black HDPE drum			
GHCL-3101-78	Guanidine HCl 6M Solution, 1135 L	Active	1135.00	L	330.00	GAL	Tote			
GHCL-3101-80	Guanidine HCl 6M Solution, 200 L	Active	200.00	L	55.00	GAL	Black HDPE drum			
GHCL-3201-59	Guanidine HCl, 50 KG	Active	50.00	KG	26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3202-01	Guanidine HCl, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
GHCL-3202-50	Guanidine HCl, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3202-93	Guanidine HCl, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-01	Guanidine HCl, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
GHCL-3220-10	Guanidine HCl, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-3220-25	Guanidine HCl, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3220-50	Guanidine HCl, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GHCL-3220-51	Guanidine HCl, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3220-92	Guanidine HCl, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GHCL-3220-93	Guanidine HCl, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-94	Guanidine HCl, 2.5 KG	Active	2.50	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-3221-25	Guanidine HCl, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3224-25	Guanidine HCl, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3224-93	Guanidine HCl, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
GHCL-3225-01	Guanidine HCl, 1KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
GHCL-3225-25	Guanidine HCl, 25KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3225-92	Guanidine HCl, 100G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GHCL-3250-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-3251-59	Guanidine HCl, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3253-25	Guanidine HCl, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3253-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-3255-25	Guanidine HCl, 25KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4223-01	Guanidine Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-4223-05	Guanidine Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-10	Guanidine Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-25	Guanidine Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4223-50	Guanidine Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
GHCL-4223-92	Guanidine Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
GHCL-4223-93	Guanidine Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-5220-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-5220-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-3220-14-76	Guanidine Thiocyanate, 14.76 KG	Active	14.76 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE/ 25L EZI-I	
GTHI-3220-15	Guanidine Thiocyanate, 14.76 KG	In development	14.76 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
GTHI-3220-95	Guanidine Thiocyanate, 10 G	In development	10.00 G		60.00 ML					
GTHI-4220-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-4221-01	Guanidine Thiocyanate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GTHI-4221-05	Guanidine Thiocyanate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-10	Guanidine Thiocyanate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-25	Guanidine Thiocyanate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-4221-50	Guanidine Thiocyanate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
GTHI-4221-92	Guanidine Thiocyanate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
GTHI-4221-93	Guanidine Thiocyanate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GTHI-5203-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-5221-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
HEHH-5201-25	HEPES Hemi Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-3220-06	HEPES, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
HEPE-3220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-20	HEPES, 20 KG	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
HEPE-3220-22	HEPES, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3220-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3220-91	HEPES, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
HEPE-3220-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3220-93	HEPES, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
HEPE-3221-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-05	HEPES, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3251-01	HEPES, 1kg	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-3251-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3251-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3251-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3254-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3254-50	HEPES, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3320-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3320-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3351-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3351-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3354-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3354-50	HEPES, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-01	HEPES Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-4250-05	HEPES Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-10	HEPES Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-25	HEPES Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-4250-50	HEPES Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-92	HEPES Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-4250-93	HEPES Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
HEPE-5201-26	HEPES, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
HEPE-5203-40	HEPES, 40 KG	Active	40.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-5220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-5220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPN-5201-26	HEPES Sodium, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
HEPN-5203-50	HEPES Sodium, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPN-5220-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPN-5251-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
IHCL-4101-02	6N HCl in IPA Solution, 500ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
IHCL-4101-97	6N HCl in IPA Solution, 180 KG	Active	180.00 KG		55.00 GAL		Black HDPE drum			
KBRO-2220-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2220-05	Potassium Bromide, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
KBRO-2220-10	Potassium Bromide, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
KBRO-2220-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
KBRO-2301-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2301-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
KBRO-2302-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2302-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
LARH-6250-01	Arginine (L) Hydrochloride Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LARH-6250-05	Arginine (L) Hydrochloride Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6250-10	Arginine (L) Hydrochloride Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
LARH-6250-25	Arginine (L) Hydrochloride Biotech, 25 KG	In development	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6250-50	Arginine (L) Hydrochloride Biotech, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LARH-6250-92	Arginine (L) Hydrochloride Biotech, 100 G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LARH-6250-93	Arginine (L) Hydrochloride Biotech, 500 G	In development	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LARH-6350-01	Arginine (L) Hydrochloride Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LARH-6350-05	Arginine (L) Hydrochloride Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-10	Arginine (L) Hydrochloride Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-25	Arginine (L) Hydrochloride Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6350-50	Arginine (L) Hydrochloride Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LARH-6350-92	Arginine (L) Hydrochloride Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LARH-6350-93	Arginine (L) Hydrochloride Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LCYS-4250-10	L-Cystine Dihydrochloride, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-4250-25	L-Cystine Dihydrochloride, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-4350-25	L-Cystine Dihydrochloride, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6250-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LCYS-6250-05	Cystine (L) Hydrochloride Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6250-50	Cystine (L) Hydrochloride Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6250-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LCYS-6250-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LCYS-6350-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LCYS-6350-05	Cystine (L) Hydrochloride Biotech, 5 KG	In development	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6350-50	Cystine (L) Hydrochloride Biotech, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6350-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LCYS-6350-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LGLM-6201-01	Glutamine (L) Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LGLM-6201-05	Glutamine (L) Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-10	Glutamine (L) Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-25	Glutamine (L) Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6201-50	Glutamine (L) Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6201-92	Glutamine (L) Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LGLM-6201-93	Glutamine (L) Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LGLM-6301-01	Glutamine (L) Biotech, 1 KG	In development	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LGLM-6301-05	Glutamine (L) Biotech, 5 KG	In development	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6301-10	Glutamine (L) Biotech, 10 KG	In development	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6301-25	Glutamine (L) Biotech, 25 KG	In development	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6301-50	Glutamine (L) Biotech, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6301-92	Glutamine (L) Biotech, 100 G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LGLM-6301-93	Glutamine (L) Biotech, 500 G	In development	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LHMM-6250-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LHMM-6250-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6250-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
LHMM-6250-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LHMM-6250-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6250-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LHMM-6250-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LHMM-6350-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LHMM-6350-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LHMM-6350-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6350-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LHMM-6350-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESH-3201-50	MES Hydrate, 50KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESH-3201-92	MES Hydrate, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESH-3250-01	MES Hydrate, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESH-3250-25	MES Hydrate, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-3350-01	MES Hydrate, 1KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESH-4250-01	MES Hydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESH-4250-05	MES Hydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-10	MES Hydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-25	MES Hydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-4250-50	MES Hydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESH-4250-92	MES Hydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESH-4250-93	MES Hydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESH-9250-91	MES Hydrate, 50 G	In development	50.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3220-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESM-3220-05	MES Monohydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-12	MES, Monohydrate, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3220-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESM-3220-91	MES Monohydrate, 50G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3220-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3220-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3221-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3221-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3222-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESM-3222-03	MES Monohydrate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36		
MESM-3222-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3222-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESM-3222-91	MES Monohydrate, 50G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3222-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3222-93	MES, Monohydrate, 500G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3223-05	MES, Monohydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-06	MES, Monohydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
MESM-3250-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
MESM-3250-12	MES, Monohydrate, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3250-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESM-3250-91	MES, Monohydrate, 250G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3250-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3250-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3250-95	MES, Monohydrate, 10G	In development	10.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3251-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3251-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3320-25	MES, Monohydrate, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3320-92	MES, Monohydrate, 100G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-3321-10	MES Monohydrate, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3321-93	MES, Monohydrate, 500G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-3350-25	MES, Monohydrate, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3350-92	MES, Monohydrate, 100G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-4220-01	MES, Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESM-4220-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4220-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4220-26	MES, Monohydrate, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
MESM-4221-01	MES, Monohydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MESM-4221-05	MES, Monohydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-10	MES, Monohydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-25	MES, Monohydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4221-50	MES, Monohydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESM-4221-92	MES, Monohydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-4221-93	MES, Monohydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESM-5201-26	MES, Monohydrate, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
MESM-5202-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-5220-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-5220-92	MES, Monohydrate, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MESM-5220-93	MES, Monohydrate, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MESN-5220-10	MES Sodium, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MESN-5220-26	MES Sodium, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
MESN-5220-50	MES Sodium, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MESN-5220-93	MES Sodium, 300 G	Active	300.00 G		1000.00 ML		White HDPE bottle, blue lid			
MOPS-3201-10	MOPS, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-3201-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-01	MOPS, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MOPS-3220-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-3220-50	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MOPS-3220-93	MOPS, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MOPS-4220-05	MOPS, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-01	MOPS Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
MOPS-4221-05	MOPS Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-10	MOPS Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-25	MOPS Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-4221-50	MOPS Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-4221-92	MOPS Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MOPS-4221-93	MOPS Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MOPS-5201-26	MOPS, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
MOPS-5202-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-5203-10	MOPS, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-5220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-5220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE round pail, blue lid			
MOPS-5221-92	MOPS, 100 G	In development	100.00 G		250.00 ML		White HDPE round pail, blue lid			
NACL-3150-02	Sodium Chloride 5M Solution, 1 L	Active	1.00 L		1.00 L		Tepoflex sample			
NACL-3150-69	Sodium Chloride 5M Solution, 200 L	Active	200.00 L		58.00 GAL		Blue HDPE drum with metal ring for BPC		BPC bag	
NACL-3150-76	Sodium Chloride 5M Solution, 1025 L	Active	1025.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-3150-76	Sodium Hydroxide 0.5N Solution, 1005 L	Active	1005.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-3153-75	Sodium Hydroxide 1N Solution, 1005 L	In development	1005.00 L		275.00 GAL		Tote			
NAHY-3153-76	Sodium Hydroxide 1N Solution, 1005 L	Active	1005.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-4122-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4122-72	Sodium Hydroxide 10N Solution, 200 L	Active	200.00 L		55.00 GAL		Natural HDPE drum			
NAHY-4122-75	Sodium Hydroxide 10N Solution, 940 L	Active	940.00 L		275.00 GAL		Tote			
NAHY-4122-96	Sodium Hydroxide 10N Solution, Bulk	Active	13059.00 L		13059.00 L		Bulk tanker			
NAHY-4123-70	Sodium Hydroxide 2N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4123-71	Sodium Hydroxide 2N Solution, 1 L X 6 combo pack	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4123-73	Sodium Hydroxide 2N Solution, 4 x 3.78 L	In development	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4124-78	Sodium Hydroxide 5N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4126-71	Sodium Hydroxide 10N Solution Biotech, 6 X 1 L	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4126-72	Sodium Hydroxide 10N Solution Biotech, 200 L	Active	200.00 L		55.00 GAL		Black HDPE drum			
NAHY-4126-73	Sodium Hydroxide 10N Solution Biotech, 4 X 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4126-75	Sodium Hydroxide 10N Solution Biotech, 1000 L	Active	1000.00 L		275.00 GAL		Tote			
NAHY-4126-78	Sodium Hydroxide 10N Solution Biotech, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4150-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4150-71	Sodium Hydroxide 10N Solution, 1 L X 6	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4150-73	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4150-78	Sodium Hydroxide 10N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4150-79	Sodium Hydroxide 10N Solution, 942.5 L	Active	942.50 L		275.00 GAL		Tote			
NAHY-4150-81	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4150-83	Sodium Hydroxide 10N Solution, 4 x 4 L	In development	16.00 L		16.00 L		4 x 4L combo pack			
NAHY-4151-70	Sodium Hydroxide 8N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NDEC-3201-01	Sodium Decanoate, 1kg	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
NDEC-3201-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
NDEC-3201-92	Sodium Decanoate, 100 gr	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
NDEC-3201-94	Sodium Decanoate, 250 gr	Active	250.00 G		1500.00 ML		White HDPE bottle, blue lid			

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
NDEC-3220-01	Sodium Decanoate, 1kg	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
NDEC-3220-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3220-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3220-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
NDEC-3220-91	Sodium Decanoate, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
NDEC-3220-92	Sodium Decanoate, 100 G	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
NDEC-3320-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-4201-01	Sodium Decanoate Biotech, 1 KG	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
NDEC-4201-05	Sodium Decanoate Biotech, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-4201-10	Sodium Decanoate Biotech, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
NDEC-4201-27	Sodium Decanoate Biotech, 25 KG	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
NDEC-4201-92	Sodium Decanoate Biotech, 100 G	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
NDEC-4201-93	Sodium Decanoate Biotech, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
THCL-3203-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3203-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3203-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3203-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3220-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3220-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3220-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3220-50	Tris Hydrochloride, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3220-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
THCL-3220-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3220-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3221-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3221-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3250-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-28	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3250-41	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-42	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
THCL-3250-91	Tris Hydrochloride, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3253-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
THCL-3253-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3253-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3253-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3253-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3254-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3254-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3254-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3254-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3254-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3254-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3254-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3257-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3257-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3259-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3259-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3259-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
THCL-3259-91	Tris Hydrochloride, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3259-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3259-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3260-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3260-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-4251-01	Tris Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-4251-05	Tris Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-10	Tris Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-25	Tris Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-4251-50	Tris Hydrochloride Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-4251-92	Tris Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-4251-93	Tris Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRED-3250-03	Trehalose Dihydrate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-05	Trehalose Dihydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3250-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3250-50	Trehalose Dihydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
TRED-3252-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3252-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3252-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3253-06	Trehalose Dihydrate, 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3253-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-05	Trehalose Dihydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-06	Trehalose Dihydrate Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-4250-10	Trehalose Dihydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRED-4250-25	Trehalose Dihydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-4250-50	Trehalose Dihydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
TRIS-1201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-1201-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-1201-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1201-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-1201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-1220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1220-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-1220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-2255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-2255-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-2255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2257-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-2257-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-2257-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2257-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-2257-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3220-50	Tris/Tromethamine, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3251-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3251-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3251-95	Tris/Tromethamine, 80 G	Active	80.00 G		1000.00 ML		White HDPE bottle, blue lid	9 x 15	LDPE	
TRIS-3252-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3252-60	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (2x) 3	LDPE/(2x) Tyve	(2x) 5 x 8 Unit
TRIS-3254-01	Tris/Tromethamine, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3254-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3254-50	Tris/Tromethamine, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3254-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRIS-3254-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3254-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3255-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-07	Tris/Tromethamine, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-12	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3255-41	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-42	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3255-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
TRIS-3255-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-91	Tris/Tromethamine, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-94	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3257-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-01	Tris (Tromethamine) Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-4223-05	Tris (Tromethamine) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-10	Tris (Tromethamine) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-25	Tris (Tromethamine) Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-4223-50	Tris (Tromethamine) Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-4223-92	Tris (Tromethamine) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-4223-93	Tris (Tromethamine) Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-5201-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5203-59	Tris/Tromethamine, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5204-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
TRIS-5207-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-5222-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-5223-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
URAC-4201-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URAC-4202-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-91	Uracil, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
URAC-4202-92	Uracil, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URAC-4202-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			
URAC-4250-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URAC-4250-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			
URAC-4301-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
URAC-4301-06	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(5x) 24 x 36	LDPE	
UREA-3120-02	Urea 6M Solution, 500 ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
UREA-3120-76	Urea 6M Solution, 750 L	In development	750.00 L		1000.00 L		BPC shell	1000L	BPC bag	
UREA-3120-78	Urea 6M Solution, 1200 L	Active	1200.00 L		330.00 GAL		Tote		LDPE	

Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
UREA-3220-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-3220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3220-50	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
UREA-3220-89	Urea, 200g	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-3220-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-3220-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-3222-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-3222-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3223-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3223-58	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3	LDPE/Tyvek	5 x 8 Unit
UREA-3250-01	Urea, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
UREA-3250-12	Urea, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	4 Unit
UREA-3250-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3250-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-3250-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-4221-01	Urea Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
UREA-4221-05	Urea Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-10	Urea Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-25	Urea Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
UREA-4221-50	Urea Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
UREA-4221-92	Urea Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-4221-93	Urea Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-5203-05	Urea, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-5203-15	Urea, 15 KG	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-5220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-3250-01	Uridine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URID-3250-10	Uridine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-3250-25	Uridine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-3250-91	Uridine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-3250-92	Uridine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URID-3250-93	Uridine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-4250-01	Uridine Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URID-4250-05	Uridine Biotech, 5KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-10	Uridine Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-25	Uridine Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-4250-50	Uridine Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
URID-4250-92	Uridine Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URID-4250-93	Uridine Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-4250-95	Uridine Biotech, 50 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
WAFI-3150-68	Water for Injection (WFI), 200L	In development	200.00 L		200.00 L		BPC shell	200L	Typoflex BPC b	
WAFI-3150-74	Water for Injection (WFI), 3.78L	Active	3.78 L		3.78 L		Non-fluorinated HDPE bottle			