

2/16/2026

Version 1.2

3.5 GAL / 5KG PAIL SPECIFICATIONS



3.5 Gallon / 5kg HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with hole(s) for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	11.0" \pm 0.5"
Inside Top Diameter	11.25" \pm 0.5"
Outside Height, with Lid	12.5" \pm 0.5"
Outside Bottom Diameter	10.5" \pm 0.5"
UN Rating	1H2/Y19/S/YR

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.